


BLADDER CANCER MANAGEMENT

Bipolar Resection and HD-NBI



NMIBC – DIAGNOSIS AND TREATMENT


The Challenge in NMIBC Management

Reviewing the literature, 50–70% of NMIBC patients suffer from recurrences after the initial treatment in the form of superficial tumors¹. Missing small papillary tumors or CIS is an acknowledged risk and occurs at an estimated rate of 10–20%². Superficial tumors of the mucosa are typically indicated by specific capillary and vessel structures. Identifying these structures – especially small tumors under white-light conditions – is possible but could be improved.


Where HD-NBI Is available:

- Rigid and flexible cystoscopy
- TUR-B resection
- Flexible ureteroscopy for upper tract tumor management

Images from a CYF-VH Flexible Cystoscope in HD and HD-NBI



Excellent view of mucosal and vessel structures even with HDTV white light in flexible video cystoscopy





HD-NBI filters the white light and improves the visual contrast of mucosal and vessel structures

¹ Sylvester et al., *Intravesical bacillus Calmette-Guerin reduces the risk of progression in patients with superficial bladder cancer: a meta-analysis of the published results of randomized clinical trials*, *J Urol* 2002; 168:1964–1970

² Jichlinski et al., *Fluorescence cystoscopy in the management of bladder cancer: a help for the urologist!*, *Urol. Int.* 2005; 74: 97–101



What Is NBI?

NBI is an optical image-enhancement technology that improves the visibility of vessels and other tissues on the mucosal surface. Narrow-band illumination, which is strongly absorbed by hemoglobin and penetrates only the surface of tissues, is ideal for enhancing the contrast between the two structures. As a result, under narrow-band illumination, capillaries on the mucosal surface are displayed in brown on the monitor, and veins in the submucosa are displayed in cyan.



What Is PDD?

When a photosensitive marker has been introduced transurethrally into the bladder, the inner surface of the bladder absorbs the drug over a period of 90–120 minutes and converts it into an endogenous pigment. This pigment is then selectively deposited in a tumor and, under blue excitation light, will emit red fluorescence. Nevertheless, in this condition, a sufficient contrast of the red fluorescence against the blue background is hard to obtain because the red fluorescence is too weak compared with the blue light. To emphasize this fluorescence, a yellow filter exclusively designed for PDD is built into the scope. As a result, the red fluorescence can be observed with excellent contrast.




CYF-VH: THE WORLD'S FIRST FLEXIBLE HD-NBI CYSTO-NEPHRO VIDEOSCOPE



HD and Full-Screen-Capabilities – for the First Time Ever in Flexible Cystoscopy

- Increased light output (20% brighter) for advanced Narrow Band Imaging (NBI)
- Enhanced visibility of superficial capillary and mucosal vessel system
- Potentially higher and improved detection rate of carcinoma in situ



A 20% brighter image improves the detection and diagnosis of suspicious tissue areas

HD-NBI: The Advanced Diagnostic Option in Bladder Cancer Follow-Up

- The combination of HD resolution and brighter NBI opens a new chapter in advanced bladder cancer follow-up diagnosis using flexible cystoscopy
- Higher tumor detection rate for HD-NBI (94.7%) vs. white light (79.2%)¹

¹ Cauberg et al., Narrow Band Imaging Cystoscopy Improves the Detection of Non-Muscle-Invasive Bladder Cancer, Urology 2010; 76:658–663

HD-NBI FOR ACCURATE TUMOR RESECTION

Crystal Clear, Brilliant, and True to Life


- Details and colors are sharp and true to life
- Significantly improved color reproduction compared with previous generation
- Advanced option for HD resection of BPH and bladder cancer due to improved visibility of vessels on the mucosal surface

HD Camera Head with NBI – a Unique Combination for Bladder Cancer Resection


- Combination of HD and NBI delivers 20% brighter light compared to prior camera-head systems
- HD-NBI reduces the recurrence risk of NMIBC by at least 10% at one year¹

Rigid Cystoscopy

- Durable all-stainless-steel design with maintenance-free stopcocks
- Gentle insertion due to smooth beak design
- Engineered and manufactured in Germany



White light



Visibility of vessels and mucosal surface improved by NBI

¹ Naselli et al., A Randomized Prospective Trial to Assess the Impact of Transurethral Resection in Narrow Band Imaging Modality on Non-Muscle-Invasive Bladder Cancer Recurrence, *European Urology* 2012; 61:908–913


TURis 2.0 BIPOLAR RESECTION

The Benefits of TURis 2.0




Safety¹

- Reduced risk of TUR syndrome²
- Minimized stimulation of obturator nerve
- Extended operation times – way beyond monopolar
- Improved teaching options
- No need to reoperate




Bloodless

- Reduced perioperative blood loss due to safe bipolar hemostasis
- Plasma vaporization with the PlasmaButton provides continuous hemostasis³



Time Saving⁴

- Self-cleaning effect of loop wire through plasma activation
- Faster postoperative recovery
- Potentially decreased catheterization times



¹ Puppo P, et al. Bipolar transurethral resection in saline (TURis): outcome and complication rates after the first 1000 cases. *J Endourol.* 2009 Jul; 23(7):1145–1149

² Michielsen DP, et al. Bipolar transurethral resection in saline – an alternative surgical treatment for bladder outlet obstruction? *J Urol.* 2007 Nov; 178(5):2035–2039

³ Geavlete B, et al. Continuous vs conventional bipolar plasma vaporisation of the prostate and standard monopolar resection: a prospective, randomised comparison of a new technological advance. *Volume 113, Issue 2, pages 288–295, February 2014*

⁴ Fagerström T, et al. Complications and clinical outcome 18 months after bipolar and monopolar transurethral resection of the prostate. *J Endourol.* 2011 Jun; 25(6):1043–1049

NMIBC MANAGEMENT PLATFORM

VISERA ELITE Imaging Platform

- Patented HD-NBI filter
- Increase in brightness by 20% compared with the previous generation for optimized imaging

ESG-400 – Intelligent HF Technologies

The award-winning OES Pro resectoscopes and the established range of TURis electrodes can now be powered by the latest-generation universal HF generator – the ESG-400.

- Automated saline detection to ensure a safe procedure
- Leakage-protection sensor to permanently ensure the highest degree of safety for the user and patient
- High Power Cut Support (HPCS) – optimizing resection in saline
- Fast Spark Monitor (FSM) – constant cutting quality
- Large, illuminated touch-screen panel
- Three effect options – control the coagulation zone



VISERA ELITE PORTFOLIO OVERVIEW

Rigid Cystoscope

A20919A 17 Fr. compact cystoscope



A22003A Therapeutic cystoscope
A20972A
A20934A
O0122



A3764A 9.5 Fr. optical urethrotome
A37014A
A37028A
A37009A



A22002A Rigid optical forceps
A20713A
A20913A



A22001A SurgMaster – TURis for
WA22366A resection in saline
A22041A
A22021A



TURis Electrodes

WA22557C PlasmaButton, 12° and 30° for
plasma vaporization



WA22302D Loop, 12°, medium



WA22558C Angled loop, 12° and 30° for
TUEB (transurethral enucleation)



Flexible Videoscopes

N3828460 CYF-VH – cysto-nephro videoscope

N3828560 CYF-VHA – cysto-nephro videoscope with aspiration

N3828660 CYF-VHR – cysto-nephro videoscope
with backward angulation



Camera Head

N4488560 CH-S190-08-LB – HD urology
ultralight camera head



N2487040 OTV-S7ProH-FD – PDD
camera head



Video Processor and Light Source

N3643860 VISERA ELITE – OTV-S190
video processor



N3643950 VISERA ELITE – CLV-S190
light source



Video Processor

WB91051W HF unit ESG-400



Monitor

N3629160 OEV261H – Olympus
high-definition LCD monitor, 26"



Insufflator

N3829670 UHI-4 – insufflator



Full-HD Recording Platform

N3808460 IMH-10 – Blu-ray and HDD
recording platform



N3808560 IMH-20 – Blu-ray and HDD
recording platform



Trolley

K10021613 VISERA ELITE – WM-NP2
mobile workstation



Specifications, design, and accessories are subject to change without any notice or obligation on the part of the manufacturer.

OLYMPUS

OLYMPUS EUROPA SE & CO. KG

Postbox 10 49 08, 20034 Hamburg, Germany
Wendenstrasse 14–18, 20097 Hamburg, Germany
Phone: +49 40 23773-0, Fax: +49 40 237765
www.olympus-europa.com